




FAMED

Żywiec

Member of REINSBERG GROUP



ACCESSORIES - OPERATING THEATRE
FURNITURE



WG-29
MAYO INSTRUMENT TABLE

Made of stainless steel, the countertop with an edge prevents tools from falling. Can be rotated horizontally by 360° (with rotation lock). Hydraulically lifted by a foot lever (range: 960 - 1330 mm). Mobile base with 3 steel-encased wheels, all equipped with a brake.

Detailed information

Dimensions: 738 x 488 mm

Height: 950 - 1350 mm

Usable surface of the table top: 680 x 430 mm

Castors diameter: 75 mm



SM-01.2
INSTRUMENT TROLLEY

Made of stainless steel, equipped with a drawer, a tiltable bowl, two accessory rails, a frame for a waste bag (with a hinged cover), a glove container (enclosed on 3 sides) and a basket for catheters. The trolley has 4 wheels with bumpers, including two with a blockade. Two push handles protrude from the table top.

Detailed information

Dimensions: 800 x 480 x 830 mm

Usable surface of the table top: 380 x 430 mm

Usable surface of the drawer: 405 x 385 x 70 mm

Castors diameter: 100 mm



SM-01.3
INSTRUMENT TROLLEY

Made of stainless steel, equipped with 3 drawers, a worktop, a shelf under the drawers, a foldable top and two accessory rails. The trolley has 4 wheels with bumpers, including two with a blockade. A push handle protrudes from the top.

Detailed information

Dimensions with folded top:
530 x 570 x 900

Dimensions with unfolded top:
860 x 530 x 900

Usable surface of the table top: 430 x 385 mm

Usable surface of the drawer: 430 x 385 x 100 mm

Castors diameter: 100 mm



SM-01.4
INSTRUMENT TROLLEY

Made of stainless steel, equipped with two tops and two accessory rails. The trolley has 4 wheels with bumpers, including two with a blockade. A push handle protrudes from the top.

Detailed information

Dimensions: 1000 x 500 x 900 mm

Usable surface of the table top: 950 x 450 mm

Castors diameter: 100 mm



SM-01.5
INSTRUMENT TROLLEY

Made of stainless steel, equipped with a worktop (with a raised edge), two baskets, a frame for a waste bag (with a plastic cover) and a handle for transport (diameter 20 mm). The trolley has 4 wheels, two of them with a blockade.

Detailed information

Dimensions: 910 x 380 x 880 mm

Usable surface of the table top: 580 x 310 x 10 mm

Dimensions of the basket: 655 x 410 x 90 mm

Castors diameter: 75 mm



SM-01.6
INSTRUMENT TROLLEY

Made of stainless steel, equipped with two worktops (with a raised edge), a basket, a frame for a waste bag (with a plastic cover), a drawer and a handle for transport (diameter 20 mm). The trolley has 4 wheels, two of them with a blockade.

Detailed information

Dimensions: 880 x 400 x 880 mm

Usable surface of the upper table top: 580 x 349 mm

Usable surface of the lower table top: 590 x 330 mm

Dimensions of the basket: 360 x 595 mm

Usable surface of the drawer: 555 x 340 x 90 mm

Castors diameter: 75 mm



SM-01.7
INSTRUMENT TROLLEY

Made of stainless steel, equipped with a worktop (with a raised edge), two plastic containers, two stainless steel bowls (that can be tilted out of the trolley), a frame for a waste bag (with a plastic cover) and a handle for transport (diameter 20 mm). The trolley has 4 wheels, two of them with a blockade.

Detailed information

Dimensions: 859 x 425 x 767 mm

Usable surface of the table top: 685 x 425 x 25 mm

Bowl diameter: 220 mm

Capacity: 2,5 liters

Storage container dimensions: 400 x 300 x 75 mm

Castors diameter: 75 mm



SM-01.8
INSTRUMENT TROLLEY

Made of stainless steel, equipped with two tops (with removable trays) and a transport handle. The trolley has 4 wheels, two of them with a blockade.

Detailed information

Dimensions: 1080 x 500 x 900 mm

Usable surface of the table top: 950 x 450 x 13 mm

Castors diameter: 100 mm



WG-02.5
TREATMENT TROLLEY

Made of stainless steel, equipped with two flat tops (the bottom one with crossbars on three sides), two stainless steel bowls (they can be folded out of the trolley) and two transport handles. The trolley has 4 wheels, two of them with a blockade.

Detailed information

Dimensions: 740 x 440 x 810 mm

Usable surface of the table top: 740 x 440 mm

Bowl diameter: 220 mm

Bowl capacity: 2,5 liters

Castors diameter: 50 mm



WG-01.31
PROCEDURE TROLLEY

Made of stainless steel, equipped with a worktop (with a raised edge) and four drawers (130 mm fronts). The five-star base is equipped with double wheels in a plastic housing, including two with a blockade.


Detailed information

Dimensions: 450 x 450 x 820 mm

Usable surface of the table top: 450 x 450 mm

Base diameter: 600 mm

Castors diameter: 75 mm





WG-07.2
MEDICAL SWIVEL STOOL

Made of stainless steel, equipped with an upholstered seat (washable, resistant to disinfection). Five-star base, wheels with a steel casing (including two with a blockade), with a leg support. Height adjustment (pneumatic actuator) operated by a lever under the seat.

Detailed information

Base diameter: 480 mm

Seat diameter: 380 mm

Height (version W1): 590 - 890 mm

Height (version W2): 510 - 710 mm

Castors diameter: 50 mm



WG-07.5
MEDICAL SWIVEL STOOL WITH BACKREST

Made of stainless steel, equipped with an upholstered seat and backrest (washable, resistant to disinfection). Five-star base, wheels with a steel casing (including two with a blockade), with a leg support. Height adjustment (pneumatic actuator) operated by a lever under the seat.

Detailed information

Base diameter: 600 mm

Seat diameter: 380 mm

Height (version W1): 590 - 890 mm

Height (version W2): 510 - 710 mm

Castors diameter: 50 mm



WK-02
MOBILE INFUSION STAND

Made of stainless steel, equipped with 2 drip hooks (optional, it can be equipped with 4 hooks) and a stainless steel base. Two wheels with a lock. Stainless steel pipe column.

Detailed information

Height: 1270 - 2320 mm

Base diameter: 530 mm



Famed Żywiec Sp. z o.o.
Fabryczna 1
34-300 Żywiec, Poland

Office:
tel.: +48 33 866 62 00
fax: +48 33 475 58 90

Export:
tel.: +48 33 866 63 75
fax: +48 33 475 58 90

export@famed.com.pl
www.famed.com.pl



The products meet the requirements of the European Directive MDD 93/42/EEC concerning medical devices and the Act on medical devices.

Version 196/2024/03/2

Subject to change due to technical progress. Pictures may differ from the actual appearance of the accessories. When choosing them, we suggest contacting our specialists and we advise not to judge their functionalities based only on photos. Famed Żywiec Sp. z o.o. also makes accessories for special customer orders, whose appearance, work and capabilities may differ from those listed in the catalog. All accessories for operating tables manufactured by our company undergo full safety tests and quality controls.

Compatibility of accessories with previous models of Famed Żywiec tables and with tables from other manufacturers is given on request.