

Technical data

	SU-05 OPTIMA 5.20	SU-05 OPTIMA 5.15	SU-05 OPTIMA 5.9	SU-05 OPTIMA 5.12
Length of the table top and total width with rails (dimensions for 5 sections, standard leg and head sections)	2100 x 600 mm		2060 x 550 mm	
Minimum height, no mattress	550 mm	680 mm		
Maximum height, no mattress	1020 mm	1080 mm		
Maximum load (and SWL)	460 kg (300 kg)		460 kg (275 kg or 300 kg - optional)	
Leg section vertical adjustment	-90° to +25° (one-piece leg plate)		-90° to +30°/+45° (one-piece leg plate/ pair of leg plates)	
Leg section horizontal adjustment, standard leg plates	180°			
Back section adjustment	-40° to +80°			
Longitudinal shift	400 mm			
Lateral tilts	30°			
Trendelenburg / anti-Trendelenburg	40° / 40°			
Headrest adjustment	-50° to +55°			
Flex / Reflex	220° / 100°			
Kidney bench adjustment	150 mm			
Leg section operated by a remote control	●	●	-	-
Table lock by a remote control	●	●	○	○
Emergency manual unlock	●	●	-	-
Carbon fiber table top	●	●	○	○
Table top translucency	0,25 mm Al	0,25 mm Al	0,83 mm Al	0,83 mm Al
Base and column covers made of carbon fiber	●	●	-	-
Base and column covers made of IntelliProtectPlus™ stainless steel	○	○	●	●
Profiled, flat „T“ base	●	●	●	-
„H“ type base, with big wheels	-	-	-	●
AlterSafe™ alternative drive	-	-	-	○
● Foam mattresses (mFix™) [80 mm]	●	●	○	○
○ Polyurethane mattresses [40 mm]	-	-	●	●
○ Foam mattresses (Velcro® type mount) [80 mm]	○	○	○	○
○ Override pDetach™ panel	○	○	○	○
Available table tops - number of sections	5	5	4/5	4/5
Power source	120/230 V ~ 50/60Hz			
Ingress protection	IPX4			
Batteries	24 V DC, battery autonomy: about 25 cycles		24 V DC, battery autonomy: about 30 cycles	
Weight of the table	250 kg			

IN STANDARD

●

○

○

○

○

-

KEY:

Famed Żywiec Sp. z o.o.

Fabryczna 1
34-300 Żywiec
Poland

Office:

tel.: +48 33 866 63 75
fax: +48 33 475 58 90

export@famed.com.pl
www.famed.com.pl



Edycja 139/2025/01/8. Famed Żywiec reserves the right to modify the product and specification as part of technical progress. All illustrations and photos used in this material are for illustrative purposes only and may not reflect the finished product. The people shown in the photos are not medical professionals. They are models. The device presented in the catalog is intended for use in health care facilities by authorized persons after having read the manual (IFU).

